

Benefits

- Solo practitioner can perform ultrasound-guided procedures
- Complicated procedures are faster and easier
- Superior image stability
- Easy to use: just place the arm and let go; fine adjustments can be made at any time
- Enhances teaching by allowing instructors to release the probe and point out important details on the ultrasound screen

Specifications

- Holds transducer:
 - weight: up to 300 gr
 - thickness: 1/2 - 2 in
- Full articulation:
 - height (z) range: 1 - 9 ft
 - reach (x,y) range: 1 - 4 ft
 - rotation: 360 degrees
- Weight: 7 kg (15 lbs)
- Warranty: 1 year
- Base available with locking caster wheels or gliders



The ultraStand is easy to use – simply move the flexible arm into position and release.

Sensitive fine adjustments can be made at any time.

The ultraStand securely holds a wide variety of transducers.

Images remain stable during typical patient movement, such as breathing.

With the probe suspended by the arm, a sterile field is easily maintained.

**Wellan ultraStand^(tm) Tested. Trusted.**

The Wellan ultraStand probe positioning system works as a third arm that allows clinicians to perform image-guided procedures completely solo, resulting in major time savings and productivity improvements. The ultraStand frees assistants for more essential work, enabling ultrasound guidance that is truly handsfree.

In clinical testing since 2005, the ultraStand is trusted by clinicians who use it by themselves regularly to perform ultrasound-guided nerve blocks, breast biopsies, and central lines. The ultraStand offers the confidence of a stable image along with the extra time and control afforded by the efficacy of having both hands available to perform interventional procedures.

A product of collaboration between physicians of Dartmouth-Hitchcock Medical Center, engineers from Dartmouth's Thayer School of Engineering, and Wellan Medical's research and development team, the ultraStand utilizes multiple US and foreign patents-pending licensed from Dartmouth-Hitchcock Medical Center.

ultraStand^(tm) is made in the USA by Wellan Medical

Wellan Medical, Inc.
10 Water Street
Lebanon, NH 03766 USA

phone: +1 603 676 8601
info@wellanmedical.com
www.wellanmedical.com